


FR-pRo-FH: Filter holder for FR-pRo units

FR-pRo-FH is a filter holder that can be mounted on the light emission port of any FR-pRo unit.

It can accommodate any light filter (low-pass, high-pass, band-pass) and configure the emission spectrum accordingly.



FR-pRo-FH can accommodate any commercially available filter of either 1/2inch or 1inch diameter. The filter can be either low-pass, high-pass or band-pass.

Features & Apps

- Optical and Thickness measurements at selected spectral regimes
- Easily mounting on any FR-pRo tool
- Ideal for the evaluation of photoresists films
- Fully configurable by the usage of any commercially available low-pass, band-pass and high-pass filter